

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

Determination of Cannabinoids Using LC-MS/MS

Scope and Application

This method is used to quantitate Δ^9 -tetrahydrocannabinol (Δ^9 -THC), Δ^8 -tetrahydrocannabinol (Δ^8 -THC), Δ^9 -tetrahydrocannabinolic acid (THCA), cannabidiol (CBD), cannabidiolic acid (CBDA), and cannabinol (CBN) in alcoholic beverages (beer, wine, and distilled spirits) using deuterated Δ^9 -THC (Δ^9 -THC-d3) as the internal standard. The samples are first prepared by diluting (typically 1:100) with an acidic aqueous-organic solvent to an appropriate level and then analyzed using liquid chromatography with tandem mass spectrometry (LC-MS/MS).

Regulatory Tolerances:

As per Public Law 115-334-Dec. 20, 2018 Sec. 10113 Hemp Production, hemp is defined as "The plant *Cannabis sativa* L. and any part of that plant, including the seeds thereof and all derivatives, extracts, cannabinoids, isomers, acids, salts, and salts of isomers, whether growing or not, with a Δ^9 -THC concentration of not more than 0.3 percent on a dry weight basis." Per TTB's Hemp Policy (<https://www.ttb.gov/formulation/hemp-policy>) alcohol beverages containing hemp or hemp components should have no detectable levels of THC. Upon dilution of the raw hemp material, or processed extract, into a finished alcoholic beverage, the concentration of Δ^9 -THC and other major cannabinoids will be on the order of parts-per-million (ppm). Therefore, the focus of this method is to quantify the major cannabinoids in the low-ppm range and can be extended to higher concentrations by further sample dilution.

Levels and Limitations

Analyte	Method Detection Limit	Method Quantitation Limit	Instrumental Linear Range
Δ^9 -THC	0.07 ng/mL	0.21 ng/mL	0.5 – 100 ng/mL
Δ^8 -THC	0.06 ng/mL	0.17 ng/mL	0.5 – 100 ng/mL
THCA	0.09 ng/mL	0.27 ng/mL	0.5 – 100 ng/mL
CBD	0.07 ng/mL	0.20 ng/mL	0.5 – 100 ng/mL
CBDA	0.06 ng/mL	0.19 ng/mL	0.5 – 100 ng/mL
CBN	0.06 ng/mL	0.17 ng/mL	0.5 – 100 ng/mL

Note: The densities of the samples of interest are roughly 1.0 g/mL. Therefore, the units of $\mu\text{g/mL}$ (or ng/mL) and ppm (or ppb) may be considered interchangeable.

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
2 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

Equipment

Glassware and Supplies:

Graduated cylinders as needed
Class A volumetric flasks as needed
Positive-displacement auto pipettors capable of delivering 50 μL to 5000 μL
15 mL conical centrifuge tubes with screw cap closure
2 mL microcentrifuge tubes
2 mL autosampler vials with split-top caps
1 L amber solvent bottle (for aqueous mobile phase)

Instrumentation:

Sciex Exion UHPLC (or equivalent)
Sciex 5500+ QTrap triple quadrupole mass spectrometer (or equivalent)

Instrument Parameters (LC):

Column:	Phenomenex Kinetex C18 2.7 μm , 2.1 \times 100 mm
Guard:	Phenomenex SecurityGuard ULTRA UHPLC C18, 2.1 mm cartridge
Column temperature:	50 $^{\circ}\text{C}$
Injection volume:	5 μL
Wash Parameters:	500 μL
Flow rate:	450 $\mu\text{L}/\text{min}$
Mobile phases:	A: 0.1% formic acid in water B: 0.1% formic acid in methanol
Gradient:	0.0 min: 75% B 4.5 min: 75% B 7.5 min: 100% B 8.5 min: 100% B 9.0 min: 75% B 10.0 min: 75% B

Instrument Parameters (Ionization Source):

Ionization mode:	Electrospray (ESI) with pos/neg switching
Ion Spray voltage (IS):	5500 V (-4500 V)
Source temperature (TEM):	550 $^{\circ}\text{C}$
Curtain gas (CUR):	30
Gas 1 (GS1):	55
Gas 2 (GS2):	55
Collision Gas (CAD):	8 (Medium)

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
3 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

Instrument Parameters (MS/MS):

Analyte	Retention Time (min)	Precursor Ion (m/z)	Product Ions (m/z)	Collision Energy (V)	Declustering Potential (V)
CBD	2.4	315.1	193.1 ^a	29	114
			259.1 ^b	25	114
CBDA*	2.9	357.1	339.1 ^a	-30	-30
			245.1 ^b	-42	-30
CBN	3.7	311.1	223.1 ^a	29	126
			195.1 ^b	37	126
Δ^9 -THC-d3 (IS)	4.45	318.1	196.0 ^a	33	130
			123.1 ^b	45	130
Δ^9 -THC	4.6	315.1	193.1 ^a	31	131
			123.1 ^b	43	131
Δ^8 -THC	5.0	315.1	193.1 ^a	31	131
			123.1 ^b	43	131
THCA*	6.9	357.1	313.1 ^a	-34	-60
			245.1 ^b	-42	-60

*Negative ion mode. ^aPrimary transitions. ^bSecondary transitions.

Entrance potential (EP) globally set to 10 V (-10 V). Cell exit potential (CXP) globally set to 12 V (-12 V). Target cycle time set to 1 s for all (+ and -) MRM experiments.

Reagent and Sample Preparation and Handling

(Vendors and product numbers listed are for convenience. Equivalent products may be used.)

Methanol, Optima LC-MS grade, CAS No. 67-56-1

Water, Optima LC-MS grade or 18 Ω , CAS No. 7732-18-5

Formic acid, LC-MS grade, CAS No. 64-18-6

Cerilliant, 2 mL ampoules containing >1 mL @ 1 mg/mL in methanol or acetonitrile:

Δ^9 -Tetrahydrocannabinol (Δ^9 -THC, 98.1%)

CAS No. 1972-08-3, Cerilliant No. T-005

Δ^9 -Tetrahydrocannabinolic acid (THCA, 97.7%)

CAS No. 23978-85-0, Cerilliant No. T-093

Cannabidiol (CBD, 99.5%)

CAS No. 13956-29-1, Cerilliant No. C-045

Cannabidiolic acid (CBDA, 98.3%)

CAS No. 1244-58-2, Cerilliant No. C-144

Cannabinol (CBN, 99.4%)

CAS No. 521-35-7, Cerilliant No. C-046

Δ^8 -Tetrahydrocannabinol (Δ^8 -THC, 99.2%)

CAS No. 5957-75-5, Cerilliant No. T-032

Deuterated Δ^9 -THC (Δ^9 -THC-d3, 98.8%)

CAS No. 81586-39-2, Cerilliant No. T-011

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
4 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

Restek, 2 mL ampoules containing >1 mL @ 1 mg/mL in methanol or acetonitrile (2nd source):
 Δ^9 -Tetrahydrocannabinol (Δ^9 -THC, 98%)
CAS No. 1972-08-3, Restek No. 34067
Cannabidiol (CBD, 99%)
CAS No. 13956-29-1, Restek No. 34011

Preparation of Solutions:

- 1) Mobile phases (stable for 3 months at room temperature)
 - a. Mobile phase A: 0.1% formic acid in water (1 L)
 - i. To a 1 L graduated cylinder add approximately 250 mL water
 - ii. Add 1.0 mL formic acid and gently swirl to mix
 - iii. Add water to the 1 L mark
 - iv. Decant into a LC solvent bottle and swirl (amber glass recommended)
 - b. Mobile phase B: 0.1% formic acid in methanol (1 L)
 - i. To a 1 L graduated cylinder add approximately 250 mL methanol
 - ii. Add 1.0 mL formic acid and gently swirl to mix
 - iii. Add methanol to the 1 L mark
 - iv. Decant into a LC solvent bottle and swirl
- 2) Seal/needle washes (stable for 6 months at room temperature)
 - a. Seal/weak needle wash: 20% methanol (1 L)
 - i. To a 1 L graduated cylinder add 200 mL methanol
 - ii. Add water to the 1 L mark
 - iii. Decant into a LC solvent bottle and swirl
 - b. Strong needle wash: 100% methanol (1 L)
 - i. Add methanol directly into a LC solvent bottle
- 3) Extraction solvents (stable for 6 months at room temperature)
 - a. Extraction solvent 1: 1.0% formic acid in 75% methanol (100 mL)
 - i. To a 100 mL graduated cylinder add 25 mL water
 - ii. Add 1.0 mL formic acid and gently swirl to mix
 - iii. Add methanol to the 100 mL mark
 - iv. Decant into a glass solvent bottle and swirl
 - b. Extraction solvent 2: 0.1% formic acid in 75% methanol (100 mL)
 - i. To a 100 mL graduated cylinder add 25 mL water
 - ii. Add 100 μ L formic acid and gently swirl to mix
 - iii. Add methanol to the 100 mL mark
 - iv. Decant into a glass solvent bottle and swirl

Preparation of Standards:

- 1) Stock standard mixes (stable for 12 months, stored at -20 °C)
 - a. Positive mode stock standard mix containing 100 μ g/mL each of Δ^9 -THC, Δ^8 -THC, CBD, and CBN (10 mL) in methanol.
 - i. Transfer entire contents from the individual standard ampoules (Cerilliant) to 2 mL microcentrifuge tubes using disposable pipettes

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024

Implementation
Date:
09/04/2024

Page
5 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

- ii. Transfer 1 mL of each standard to a 10 mL volumetric flask using an auto pipette
 - iii. Fill to the line of the flask with methanol
 - iv. Shake and transfer to a 15 mL centrifuge tube
 - b. Negative mode stock standard mix containing 100 µg/mL each of THCA and CBDA (10 mL) in acetonitrile.
 - i. Transfer entire contents from the individual standard ampules (Cerilliant) to 2 mL microcentrifuge tubes using disposable pipettes
 - ii. Transfer 1 mL of each standard to a 10 mL volumetric flask using an auto pipette
 - iii. Fill to the line of the flask with acetonitrile
 - iv. Shake and transfer to a 15 mL centrifuge tube
 - c. Stock internal standard containing 100 µg/mL of Δ^9 -THC-d3 (10 mL) in methanol.
 - i. Transfer entire contents from the ampule (Cerilliant) to a 2 mL microcentrifuge tube using a disposable pipette
 - ii. Transfer 1 mL of the standard to a 10 mL volumetric flask using an auto pipette
 - iii. Fill to the line of the flask with methanol
 - iv. Shake and transfer to a 15 mL centrifuge tube
- 2) Intermediate stock standard mixes (stable for 12 months, stored at -20 °C)
 - a. Positive mode intermediate stock standard mix containing 20 µg/mL of Δ^9 -THC, Δ^8 -THC, CBD, and CBN (10 mL) in methanol.
 - i. Transfer 2 mL of positive mode stock standard mix (100 µg/mL) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - b. Negative mode intermediate stock standard mix containing 20 µg/mL of THCA and CBDA (10 mL) in acetonitrile.
 - i. Transfer 2 mL of negative mode stock standard mix (100 µg/mL) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - c. Positive mode intermediate stock standard mix containing 200 ng/mL of Δ^9 -THC, Δ^8 -THC, CBD, and CBN (10 mL) in methanol.
 - i. Transfer 100 µL of the positive mode intermediate stock standard mix (20 µg/mL) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - d. Negative mode intermediate stock standard mix containing 200 ng/mL of THCA and CBDA (10 mL) in acetonitrile.

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024

Implementation
Date:
09/04/2024

Page
6 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

- i. Transfer 100 μL of the negative mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - e. Intermediate stock internal standard containing 1 $\mu\text{g}/\text{mL}$ of $\Delta^9\text{-THC-d3}$ (10 mL) in methanol.
 - i. Transfer 100 μL of stock internal standard (100 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
- 3) Positive mode working standard mixes containing $\Delta^9\text{-THC}$, $\Delta^8\text{-THC}$, CBD, CBN, and $\Delta^9\text{-THC-d3}$ in methanol (stable for 6 months, stored at $-20\text{ }^\circ\text{C}$)
 - a. 10X Calibration level 6 (1 $\mu\text{g}/\text{mL}$, 10 mL)
 - i. Transfer 500 μL of the positive mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - b. 10X Calibration level 5 (500 ng/mL, 10 mL)
 - i. Transfer 250 μL of the positive mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - c. 10X Calibration level 4 (100 ng/mL, 10 mL)
 - i. Transfer 50 μL of the positive mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - d. 10X Calibration level 3 (50 ng/mL, 10 mL)
 - i. Transfer 2.5 mL of the positive mode intermediate stock standard mix (200 ng/mL) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - e. 10X Calibration level 2 (10 ng/mL, 10 mL)
 - i. Transfer 500 μL of the positive mode intermediate stock standard mix (200 ng/mL) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - f. 10X Calibration level 1 (5 ng/mL, 10 mL)
 - i. Transfer 250 μL of the positive mode intermediate stock standard mix (200 ng/mL) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with methanol

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
7 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

- iii. Shake and transfer to a 15 mL centrifuge tube
- 4) Negative mode working standard mixes containing THCA, CBDA, and Δ^9 -THC-d3 in acetonitrile (stable for 6 months, stored at -20 °C)
- a. 10X Calibration level 6 (1 $\mu\text{g}/\text{mL}$, 10 mL)
 - i. Transfer 500 μL of the negative mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - b. 10X Calibration level 5 (500 ng/mL, 10 mL)
 - i. Transfer 250 μL of the negative mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - c. 10X Calibration level 4 (100 ng/mL, 10 mL)
 - i. Transfer 50 μL of the negative mode intermediate stock standard mix (20 $\mu\text{g}/\text{mL}$) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - d. 10X Calibration level 3 (50 ng/mL, 10 mL)
 - i. Transfer 2.5 mL of the negative mode intermediate stock standard mix (200 ng/mL) and 1 mL of the intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - e. 10X Calibration level 2 (10 ng/mL, 10 mL)
 - i. Transfer 500 μL of the negative mode intermediate stock standard mix (200 ng/mL) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
 - f. 10X Calibration level 1 (5 ng/mL, 10 mL)
 - i. Transfer 250 μL of the negative mode intermediate stock standard mix (200 ng/mL) and 1 mL of intermediate stock internal standard (1 $\mu\text{g}/\text{mL}$) to a 10 mL volumetric flask using an auto pipette
 - ii. Fill to the line of the flask with acetonitrile
 - iii. Shake and transfer to a 15 mL centrifuge tube
- 5) Positive mode calibration standard mixes containing Δ^9 -THC, Δ^8 -THC, CBD, CBN, and Δ^9 -THC-d3 in mobile phase (stable for 24 hours at ambient temperature)
- a. For 1 mL at each calibration level of the positive mode working standard mixes:
 - i. Transfer 100 μL of the working standard mix to an autosampler vial
 - ii. Add 900 μL of 0.1% formic acid in 75% methanol
 - iii. Cap and shake

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024

Implementation
Date:
09/04/2024

Page
8 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

- 6) Negative mode calibration standard mixes containing THCA, CBDA, and Δ^9 -THC-d3 in mobile phase (stable for 24 hours at ambient temperature)
 - a. For 1 mL at each calibration level of the negative mode working standard mixes:
 - i. Transfer 100 μ L of the working standard mix to an autosampler vial
 - ii. Add 900 μ L of 0.1% formic acid in 75% methanol
 - iii. Cap and shake

Preparation of Laboratory Control Sample (LCS):

- 1) Transfer 4.85 mL of distilled spirit LCS (e.g. 40% tequila) to a 15 mL centrifuge tube
- 2) Add 50 μ L of the positive mode stock standard mix containing 100 μ g/mL each Δ^9 -THC, Δ^8 -THC, CBD, and CBN in methanol.
- 3) Add 50 μ L of the negative mode stock standard mix containing 100 μ g/mL each THCA and CBDA in acetonitrile.
- 4) Add 50 μ L of stock I.S. solution containing 100 μ g/mL Δ^9 -THC-d3 in methanol
- 5) Cap and shake
- 6) Transfer 100 μ L of spiked sample prepared in steps 1-3 to a 2 mL microcentrifuge tube
- 7) Add 900 μ L of 1.0% formic acid in 75% methanol (Extraction Solvent 1)
- 8) Cap and shake
- 9) Transfer 100 μ L of diluted spiked sample to an autosampler vial
- 10) Add 900 μ L of 0.1% formic acid in 75% methanol (Extraction Solvent 2)
- 11) Cap and shake
- 12) Repeat steps 5 – 10 to prepare the 2nd LCS sample

Preparation of 2nd Source Sample (QC):

- 1) Prepare the QC stock standard mix containing 100 μ g/mL each of Δ^9 -THC and CBD (Restek) in the same manner as stated above (10 mL).
- 2) Prepare the QC intermediate stock standard mix containing 20 μ g/mL each of Δ^9 -THC and CBD in the same manner as stated above (10 mL).
- 3) Prepare the QC working standard mix containing 100 ng/mL each of Δ^9 -THC and CBD in the same manner as stated above (10 mL).
- 4) Prepare the QC sample by transferring 100 μ L of QC working standard mix to an autosampler vial, adding 900 μ L of 0.1% formic acid in 75% methanol (Extraction Solvent 2), capping, and shaking.

Preparation of Samples:

- 1) Determine the density of the sample using a densitometer
- 2) Transfer 5 mL of sample using to a 15 mL centrifuge tube
- 3) Add 50 μ L of stock I.S. solution (100 μ g/mL Δ^9 -THC-d3 in methanol)
- 4) Cap and shake
- 5) Transfer 100 μ L of I.S.-spiked sample prepared in steps 1-3 to a 2 mL microcentrifuge tube

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
9 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

- 6) Add 900 μ L of 1.0% formic acid in 75% methanol (Extraction Solvent 1)
- 7) Cap and shake
- 8) Transfer 100 μ L of diluted spiked sample to an autosampler vial
- 9) Add 900 μ L of 0.1% formic acid in 75% methanol (Extraction Solvent 2)
- 10) Cap and shake

Samples should be prepared and analyzed on the same day (within 24 hours).

Procedures

- 1) Prepare fresh calibration standards from working standards as described in the Preparation of Standards section above.
- 2) Prepare the samples as described in the Preparation of Samples section above.
- 3) Prime the LC equilibrate the system to the initial instrument method conditions.
- 4) Inject blanks, calibrants, samples, and QC using the following recommended sequence template:
 - Solvent blank
 - Positive mode calibration levels 1 – 6
 - Negative mode calibration levels 1 – 6
 - Solvent blank
 - Two LCS
 - Solvent blank
 - QC check (2nd source sample)
 - Solvent blank
 - Sample(s)
 - Solvent blank
 - QC check (2nd source sample)
 - Solvent blank
- 5) Shut down system and properly store LC column.
- 6) Process sequence in SciexOS Analytics using processing method.
- 7) Report results as described in the Reporting Results section below.

Quality Control

- 1) The correlation coefficient (R^2) for all calibration curves is to be ≥ 0.995 . If R^2 is < 0.995 , prepare fresh calibration standard mixes and re-run injections.
- 2) Run 2 LCS samples for accuracy and precision. The values for accuracy and precision are to be within the prescribed limits for both Δ^9 -THC and CBD. If the values are outside the prescribed limits, samples are to be prepared fresh and re-run.
- 3) Run a 2nd source check standard (QC) at least once every 8 samples (including LCS). The QC check is to be 1.0 ± 0.15 ppm for both Δ^9 -THC and CBD. If the QC check is not within the prescribed limit, all standards/samples run prior to the QC check are to be re-run.

Courtesy Copy

SSD:TM:317

Rev. 2

Issue Date:
08/21/2024Implementation
Date:
09/04/2024Page
10 of 10

The colored ink stamp indicates this is a controlled document. Absence of color indicates this copy is not controlled and will not receive revision updates.

Sources of Uncertainty

- Pipetting errors
- Impure and/or contaminated standards or reagents
- Matrix interferences (overlapping signal with analyte)

Calculations

$$\text{Actual Conc. (ppm)} = \frac{\text{Calc. Conc. (ng/mL)}}{\text{Sample Density (g/mL)}} \times \frac{\text{Dilution Factor}}{1000}$$

Reporting Results

Report the results as follows:

Analytes are calculated in ng/mL and reported in ppm (XX.xx or X.xx).

If calc. conc. is <0.1 ng/mL, result is reported as "Not Detected".

If calc. conc. is ≥ 0.1 ng/mL and <0.5 ng/mL, result is reported as "Below Quantitation Limit".

If calc. conc. is ≥ 0.5 ng/mL, result is reported as "XX.xx or X.xx ppm".

If calc. conc. is >100 ng/mL, analyst must re-analyze with larger dilution factor.

Safety Notes

Anticipated waste volume for each sample preparation and UPLC run is approximately 15 mL consisting primarily of aqueous organic solvent (methanol and/or ethanol with water). Will possibly contain up to percent levels of some cannabinoids.

Required Training, Certification and Re-certification

- 1) Receive in-house LC-MS/MS training.
- 2) Initial certification is achieved by running 3 LCS replicates with results of precision and accuracy in agreement with the results of the validation package.
- 3) Periodically, analysts are retested for competency (e.g. every 5 years) and/or given proficiency test.

Revision History

Rev. 1 – Initial Revision

Rev. 2 – Addition of Δ^8 -THC as an analyte and Δ^9 -THC-d3 as the internal standard